

Cirona™ 6200 Deep Vein Thrombosis (DVT) Prevention Therapy System

The Cirona™ 6200 Deep Vein Thrombosis (DVT) Prevention
Therapy System is a pneumatic compression device that
noninvasively helps reduce the incidence of deep vein
thrombosis, a potentially life-threatening condition in which
a blood clot forms in a vein deep inside the body. If the blood
clot breaks off, it can travel to an artery in the lungs and block
blood flow.

The system consists of a pump and a pair of soft compression garments that wrap around the calf, calf-thigh, or foot on each leg or foot. Air compression inflates the garments during a 60-second automatically timed cycle (an approximate 12-second inflation period followed by a 48-second relaxation period) to mimic the body's natural walking functionality, enhancing circulation, and prevent early blood clotting.

The Cirona™ 6200 Series system is a prescription device intended to be used preventatively to increase venous blood flow in patients at risk of deep vein thrombosis due to the associated risk factors for thrombus formation during: trauma, critical care, general medicine, general surgery, as well as neurological, orthopedic, urologic, obstetric conditions and treatments.

*Caution: In the USA, Federal Law restricts this device to sale, by or on the order of a physician.





Key Features

- Extension tubing set and hospital bed hanger included
- Unilateral and Bilateral treatment available for foot, calf, or calf-thigh
- Comfortable, disposable single-use garments designed to be compatible with the other major DVT pumps
- Quick-release connectors allow for easy removal of the garments
- Treatment can be paused and started anytime
- 6+ hour battery life
- Detachable power cord for ease of transport

LB37.0004 Rev B 20150121

Specifications

Dimensions: 4"(W) X 7"(L) X 11"(H)

Weight: 5.1 lbs.

Inflation Cycle: 12 seconds
Rest Period: 48 seconds
Cycle Time: 60 seconds
Pressure: 40mmHg

Electrical: 100 - 240 VAC, 50/60Hz

Battery: Lithium Battery, 10.8 V, 4400mAh (rechargeable)

Fuse Rated: 250 VAC, 1.0 AMP (SLO-BLO)

Protection Against Electric Shock: Class II

Operation Mode: Continuous

Protection against Water: IPXO

Available Sleeves:

Model	Application	Circumference	Patient Sizing
D-110	Calf	Up to 17"	Standard
D-120	Calf	Up to 23"	Large
D-125 *	Calf	Up to 32"	Extra Large
D-135 *	Calf - Thigh	Up to 22"	Small
D-130	Calf - Thigh	Up to 28"	Standard
D-140	Calf - Thigh	Up to 36"	Large
D-145 *	Calf - Thigh	Up to 42"	Extra Large
F-15 *	Foot	Up to 13"	Small
F-10	Foot	Up to 16"	Standard
Extension Tubing Set	6' tube extender to connect the garments to the device.		

 $\,\star\,$ These sizes are available on the order confirmation only.



The Cirona™ 6200 Series should NOT be used in the

- following conditions:
- Severe atherosclerosis or other ischemic vascular diseases
- · Suspected or known acute deep vein thrombosis
- · Severe congestive cardiac failure
- · Existing pulmonary edema
- Existing pulmonary embolisms
- · Extreme deformity of the limbs
- Any local skin or tissue condition in which the garments would interfere, such as:
 - Gangrene
 - Untreated or infected wounds
 - Recent skin grafts
 - · Dermatitis, etc.
- . Known presence of malignancy in the legs
- Limb infections, including cellulitis, that have not received antibiotic coverage
- · Presence of lymphangiosarcoma





Hanstronics Medical Equipments Pvt Ltd. 1551 Aggarwal Metro Heights, Netaji Subhash Place, Pitampura, New Delhi-110034, India.

Mobile: +91 9811011134 Ph: +91 11 47076038, 42060250

Email: pankaj@hanstronics.com Website: www.hanstronics.com