



extriCARE® 3600 Negative Pressure Wound Therapy (NPWT) System

The **extriCARE® 3600 Negative Pressure Wound Therapy (NPWT)** device is a high capacity and high flow pump designed for use in hospitals, wound care centers, and other medical facilities.

The **extriCARE® 3600 NPWT pump** is a portable, battery powered pump which may promote wound healing through the drainage and removal of wound exudates, infectious material, and tissue debris from the wound bed using continuous and/or intermittent negative pressure.

Indications for Use

The **extriCARE® 3600 Negative Pressure Wound Therapy System** is indicated for wound management via the application of negative pressure to the wound by the removal of wound exudate, infectious materials, and tissue debris from the wound bed. The **extriCARE® 3600 Negative Pressure Wound Therapy System** is indicated for the following wound types: chronic, acute, traumatic, subacute and dehiscent wounds, partial-thickness burns, ulcers (such as diabetic or pressure), flaps and grafts.

Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

Key Features:

- Long-Lasting Battery: 24 Hours
- Rotatable Back Mount for Securing to Hospital Bed Rails, IV Poles, etc.
- LCD Screen, Intuitive User Interface, and Easy Pressure Setting (40-200mmHg)
- Continuous Mode or Intermittent Mode
- Alarms: Leakage, Canister Full, Canister Installation Error, High Voltage, Low Battery, Blockage





Device Specifications:

Weight: 2.87 lbs (1.3 kg)

Mobility: Portable

Battery Type: Lithium (rechargeable)

Vacuum Modes: Continuous or intermittent

Dimensions:

Length: 6.7 (17cm)

Depth: 4.3 (11cm)

Height: 5.1 (13cm)

Input:

100-240 V

50/60 Hz

1.5 A

Operating Conditions:

Temperature: +5° C 40° C

Humidity: 15% 93%

Pressure Options: 40mmHg 200mmHg

Charging Time: < 5 hours

Barometric Pressure Range: 800 hPa 1060 hPa

Storage Temperature: -25°C - +70°C

Altitude Range: < 2000 m

Ingress Protection: IP21

Protection Against Electric Shock: CLASS II

Patient Protection: Type BF

Precautions & Warnings

Extra precaution and care should be taken when using

Negative Pressure Wound Therapy in the following conditions:

- Bleeding/Hemorrhaging
- Hemostasis, anticoagulants, and platelet aggregation inhibitors
- Friable blood vessels or infected vascular anastomosis
- Exposed bones, vessels, ligaments or tendons
- Sharp edges (bone fragments)
- Not to be used in MRI, hyperbaric chamber or defibrillation environment
- Infected wounds and osteomyelitis
- Patient size and weight
- Spinal cord injury
- Enteric Fistulas
- Circumferential Dressing application
- Use near the vagus nerve (Bradycardia)
- Protection of the periwound skin
- Intermittent versus continuous mode

Contraindications

The use of the extriCARE® 3600 pump

is contraindicated in the presence of:A

- Exposed vessels, organs, or nerves.
- Anastomotic sites.
- Exposed arteries or veins in a wound.
- Fistulas, unexplored or non-enteric.
- Untreated osteomyelitis.
- Malignancy in the wound.
- Excess amount of necrotic tissue with eschar.
- Wounds which are too large or too deep to be accommodated by the dressing.
- Inability to be followed by a medical professional or to keep scheduled appointments.
- Allergy to urethane dressings and adhesives.
- Use of topical products which must be applied more frequently than the dressing change schedule allows.



For more information Contact

Hanstronics Medical Equipments Pvt. Ltd.

1551, Aggarwal Metro Heights, Plot No. E-5, 15th Floor, Netaji Subhash Place, Pitampura, New Delhi-34

Ph.: 011-47076038, 42060250 | Email : hmepl@hanstronics.com | Web : www.hanstronics.com