



extriCARE® Negative Pressure Wound Therapy (NPWT) System Foam Dressings

The extriCARE® NPWT system incorporates the use of extriCARE® pump models with Devon Medical's foam dressings.

These dressings utilize a foam pad under an adhesive film that can be applied directly to the wound bed. A drainage tube connects each dressing with a canister that is connected to the pump. The canister collects exudates and bodily fluid that is removed from the wound bed.

Application of the foam dressing is very quick and easy, allowing wound care professionals and caregivers to decrease application time.

Foam dressings are available in various shapes and sizes and designed for use with wounds on all parts of the body.

extriCARE® Foam Dressings	Type of Dressing	Dressing Size	Foam Pad Size
	EC2400-FOAM-S-E	Small Foam Dressing	100mm x 75 mm (3.9" x 3")
	EC2400-FOAM-L-E	Large Foam Dressing	250mm x 160mm (9.8" x 6.3")