

AIROS™ 6 SEQUENTIAL COMPRESSION DEVICE

Reliable, Intuitive 6-Chamber Device
Serving the Medicare Population



The FDA-cleared AIROS 6 Sequential Compression Device is a compression therapy device utilized for the treatment and management of lymphedema and venous insufficiencies. The device provides effective intermittent pneumatic compression (IPC) therapy for treatment of the legs and arms.

The AIROS 6 provides gradient sequential compression that increases blood flow and encourages extracellular fluid clearance. The device provides cycles of compressed air at certain adjustable pressures and then inflates and deflates garments to move fluid.

KEY FEATURES

- **Billing:** Medicare patients (Code E0651) or private insurance
- **Ease of Use Promoting Patient Compliance:** Large, clear LCD screen and one-button operation for starting and stopping therapy
- **Physical Benefits :** Reduce pain and swelling, increase blood flow, move lymphatic fluid
- **Various Sizes of Comfortable Garments:** 8 full leg sizes, 2 half leg sizes, and 2 arm sizes available to comfortably fit patients (all utilize 6 chambers)
- **Lymphedema Preparation Treatment (LPT) Modes**
Optional Clearing Mode (mimics lymphatic drainage); or Focal Mode (fibrotic tissue relief) can be used before main Gradient Mode therapy

SPECIFICATIONS

Dimensions	10.2" W x 10.5" D x 4.3" H
Weight	6.4 lbs.
Primary Operating Mode	Gradient
Secondary Operating Modes (LPT)	Clearing or Focal
Pressure Range (Gradient Mode)	30-80mmHg (increments of 1mmHg)
Cycle Time (Inflation)	40 or 60 seconds
Deflation Time	12 seconds
Therapy Time	15, 30, 45 or 60 minutes
Warranty (Device)	3 years
Alarms (Error Messages)	3
Air Blocker (Unilateral Therapy)	Included (1)
Electrical	120 VAC, 60Hz, 100VA MAX

AIROS 6-CHAMBER COMPRESSION GARMENTS

Durable, Comfortable Garments
for Effective Therapy



The AIROS™ 6-Chamber Compression Garments are intended for use with the FDA-cleared AIROS 6 Sequential Compression Device. Together, this medical technology provides compression therapy for patients living with lymphedema and venous insufficiencies.

The compression garments were designed with patient comfort and clinical effectiveness in mind. From durable material to precise sizing, our garments were carefully engineered to support patients with various needs.

KEY FEATURES

- **Various Sizes:** 8 full leg sizes, 2 half leg sizes, and 2 arm sizes available to comfortably fit patients (all utilize 6 chambers)
- **Ease of Use:** Pull handles and O-ring zippers to facilitate easy application and removal
- **High-Quality Material:** Durable, easy-to-clean material also allows for simple cleaning

AIROS™ 6-Chamber Compression Garments Sizing Chart

LEG GARMENTS					
Garment	Model	Inseam	Circumference (Ankle)	Circumference (knee)	Circumference (Thigh)
Small Leg	6-SL01	66cm (26")	53cm (21")	64cm (25")	76cm (30")
Small Wide Leg	6-SWL01	66cm (26")	61cm (24")	74cm (29")	89cm (35")
Small Extra Wide Leg	6-SXWL01	66cm (26")	71cm (28")	89cm (35")	102cm (40")
Medium Leg	6-ML01	74cm (29")	56cm (22")	69cm (27")	79cm (31")
Medium Wide Leg	6-MWL01	74cm (29")	64cm (25")	79cm (31")	89cm (35")
Med. Extra Wide Leg	6-MXWL01	74cm (29")	74cm (29")	91cm (36")	102cm (40")
Large Leg	6-LL01	91cm (36")	53cm (21")	69cm (27")	79cm (31")
Large Wide Leg	6-LWL01	91cm (36")	56cm (22")	79cm (31")	102cm (40")
ARM GARMENTS					
Medium Arm	6-MA01	74cm (29")	56cm (22")	48cm (19")	33cm (13")
Large Arm	6-LA01	81cm (32")	58cm (23")	51cm (20")	36cm (14")

www.airosmedical.com



For more information, please contact :
Hanstronics Medical Equipments Pvt. Ltd.
1551, Aggarwal Metro Heights, Plot No. E-5,
Netaji Subhash Place, Pitampura,
New Delhi-110034, India

+91 11 47076038, 42060250, 9811011134
pankaj@hanstronics.com
www.hanstronics.com